Instructions For Use

SAS-I Immunofix Cat. No. 200300

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INTENDED PURPOSE

The SAS-I IFE kit is intended for the separation and identification of monoclonal gammopathies by agarose gel electrophoresis.

Immunofixation electrophoresis (IFE) is a two stage procedure using high resolution agarose electrophoresis in the first stage, and immunoprecipitation in the second phase.

The greatest demand for IFE is in the clinical laboratory, where it is used primarily for the detection of monoclonal gammopathies. A monoclonal gammopathy is a primary disease state in which a single clone of plasma cells produce elevated levels of an immunoglobulin of a single class and type. Such immunoglobulins are referred to as monoclonal proteins, M-proteins or paraproteins. Their presence may be of a benign nature or of uncertain significance. In some cases, they are indicative of a malignancy, such as multiple myeloma or Waldenström's macroglobulinaemia. Differentiation must be made between polyclonal and monoclonal gammopathies, as polyclonal gammopathies are a secondary disease state due to clinical disorders such as chronic liver disease, collagen disorders, rheumatoid arthritis and chronic infection.

Urinary proteins are derived primarily from plasma proteins that filter through the kidney. The appearance of abnormal plasma proteins in the urine is of great value in evaluating renal function. The appropriate study of proteinuria should include quantitative and qualitative assessment of the type and amount of proteins excreted^{1,5}. The combination of electrophoretic separation of urine proteins, coupled to the identification of specific protein types by immunoprecipitation allows the differentiation of several types of proteinuria - Physiological, Glomerular (selective and non-selective), Tubular and proteinuria associated with Dysglobulinaemias^{1,5}.

Alfonso first described immunofixation in the literature in 1964⁶. Alper and Johnson published a more practical procedure in 1969, and published a number of studies utilising this technique⁷⁻⁹. Immunofixation has been used as a procedure for the investigation of immunoglobulins since 1976¹⁰⁻¹¹. The SAS-I IFE kit separates serum proteins according to charge in an agarose gel. The proteins are then incubated with monospecific antisera, washed and stained to allow visualization of the immunoprecipitate for qualitative interpretation.

WARNINGS AND PRECAUTIONS

All reagents are for in-vitro diagnostic use only. Do not ingest or pipette by mouth any kit component. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information.

COMPOSITION

1. SAS-I IFE Gel

Contains agarose in a Tris / Barbitol buffer with thiomersal and sodium azide as preservative. The gel is ready for use as packaged.

2. Acid Violet Stain Concentrate

Contains concentrated Acid Violet stain. Dilute the contents of the bottle to 700ml with purified water. Stir overnight and filter before use. Store in a tightly stoppered bottle.

3. Destain Solution Concentrate

Dilute the contents of Destain A to 1 litre with purified water. Then add the contents of Destain B and add a further 1 litre of purified water, slowly.

4. Wash Solution

Contains 100ml of concentrated Wash Solution. Dilute 20ml to 1 litre with saline solution.

5. Sample Diluent

Contains Tris / Barbitol buffer with bromophenol blue and sodium azide as preservative. The diluent is ready for use as packaged.

6. SAS-I IFE Antisera Kit

Contains SP protein fixative (containing acetic acid and sulphosalicylic acid), and monospecific antisera to human immunoglobulins - IgG, IgA, IgM, Kappa Light Chain (free & bound) and Lambda Light Chain (free & bound). All antisera contain sodium azide as a preservative. The antisera are ready for use as packaged.

7. Other Kit Components

Each kit contains Instructions For Use and sufficient Blotters B, C, D, X and Blotter Combs to complete 10 gels.

STORAGE AND SHELF-LIFE

1. SAS-I IFE Gel

Gels should be stored at 15...30°C and are stable until the expiry date indicated on the package. DO NOT REFRIGERATE OR FREEZE. Deterioration of the gel may be indicated by 1) crystalline appearance indicating the gel has been frozen, 2) cracking and peeling indicating drying of the gel or 3) visible contamination of the agarose from bacterial or fungal sources.

2. Acid Violet Stain

The stain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted stain solution is stable for 6 months at 15...30°C. It is recommended to discard used stain immediately to prevent depletion of staining capability. Poor staining performance may indicate deterioration.

3. Destain Solution

The destain concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted destain solution is stable for 6 months at 15...30°C.

4. Wash Solution

The Wash Solution concentrate should be stored at 15...30°C and is stable until the expiry date indicated on the label. Diluted wash solution is stable for 6 months at 15...30°C. Cloudiness may indicate deterioration.

5. Sample Diluent

The Sample Diluent should be stored at 15...30°C and is stable until the expiry date indicated on the label. Cloudiness may indicate deterioration.

6. SAS-I IFE Antisera Kit

The antisera kit should be stored at 2...6°C and is stable until the expiry date indicated on the label. Particulate contamination or cloudiness may indicate deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

Cat. No. 210200 Sample Applicator Blades (1 x 10)

Cat. No. 210300 Sample Applicator Blades (5 x 10)

Cat. No. 210100 Disposable sample cups (100)

Cat. No. 5014 Development Weight (2.3kg)

Cat. No. 3100 REP Prep

Drying oven with forced air capable of 60...70°C.

Saline solution (0.85% NaCl)

The following items are not required for standard serum IFE but may be required for further investigation and Urine IFE investigation:

Cat. No. 220100 Antiserum to Human Urine Total Protein (2ml)

Cat. No. 220200 Antiserum to Human Urine Micro Proteins (2ml)

Cat. No. 220300 Antiserum to Human Urine Macro Proteins (2ml)

Cat. No. 220400 Antiserum to Human GAM Proteins (2ml)

Cat. No. 220700 Antiserum to Human Free & Bound Kappa Light Chain (2ml)

Cat. No. 220800 Antiserum to Human Free & Bound Lambda Light Chain (2ml)

Cat. No. 220500 Antiserum to Human Free Kappa Light Chain (2ml)

Cat. No. 220600 Antiserum to Human Free Lambda Light Chain (2ml)

Cat. No. 220900 Antiserum to Human Urine Pentavalent/Albumin (2ml)

Cat. No. 221000 Antiserum to Human Urine Pentavalent (2ml)

Cat. No. 9249 Antiserum to Human IgD

Cat. No. 9250 Antiserum to Human IgE

Cat. No. 9400 IFE Control Kit

SAMPLE COLLECTION AND PREPARATION

Freshly collected serum is the specimen of choice. Samples can be stored at 15...30°C for up to 4 days, 2...6°C for up to 2 weeks or 6 months at -20°C⁵.

Urine samples should be applied neat. Serum samples should be diluted according to the following table using Sample Diluent:

Monoclonal Concentration	SP Lane	G, A, M, K, λ
<3g/L	1+2	Neat
3-10g/L	1+2	1+4
10-25g/L	1+2	1+9
25-50g/L	1+2	1+19

STEP BY STEP PROCEDURE

- Pipette 35µl of the sample into the appropriate well of the SAS-I sample tray or disposable sample cups.
 - SAS-I & SAS-I Plus users:** Carefully place the sample tray onto the applicator drawer. Ensure that the tray is pushed firmly down into position.
 - SAS-3 users:** Carefully locate the sample tray using the sample base locating pins. Ensure that the tray is positioned securely.
- Remove the gel from the packaging and:
 - SAS-I users:** place the gel in the SAS-I, agarose side up, aligning the positive and negative sides with the corresponding electrode posts.

- ii) **SAS-1 Plus users:** dispense 400 μ L of REP Prep onto the heat sink. Place the gel onto the heat sink, agarose side up, aligning the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
- iii) **SAS-3 users:** place the alignment guide onto the pins and dispense 400 μ L of REP Prep onto the centre of the chamber. Place the gel into the chamber agarose side up, using the guide, align the positive and negative sides with the corresponding electrode posts, taking care to avoid air bubbles under the gel.
3. Blot the surface of the gel with a blotter C, discard the blotter.
4. i) **SAS-1 users:** attach the electrodes onto the top side of the electrode posts so that they are in contact with the gel blocks.
 ii) **SAS-1 Plus users:** (as above). Place the cover over the gel and electrodes and press firmly for 5 seconds to ensure contact.
 iii) **SAS-3 users:** attach the electrodes onto the the electrode posts so that they are in contact with the gel blocks.
5. Place 2 applicator blade assemblies in position on the instrument, (**SAS-3 users: slot A and 10**).
6. Perform the Immunofix electrophoresis:
 - i) **SAS-1 users:** 80 volts, 20 mins
 - ii) **SAS-1 Plus users:** Electrophoresis: 100 volts, 18 mins, 20°C
 Incubation Step 1: 8 mins, 37°C (incubate)
 Incubation Step 2: 8 mins, 40°C (D blot)
 - iii) **SAS-3 users:**

Step	Time (mm:ss)	Temperature (°C)	Voltage	Other
Load Sample	00:30	21		Speed 1
Apply Sample	00:30	21		Speed 1*
Electrophoresis	17:00	21	100	
Apply antisera	10:00	21		
Insert Combs	02:00	21		
Blotter D	05:00	40		
Dry	08:00	54		

* Use Location 2

NOTE 1: For Serum immunofixation, 1 sample application is required. For Urine immunofixation, 10 sample applications are required. Remove the gel blocks prior to drying.

NOTE 2: Urine and serum samples can be run simultaneously on one gel. However, each individual row must contain one type of sample only (eg. top row = urine, bottom row = serum). If a combination of serum and urine are to be used on the same gel, place both blades onto the SAS-1plus, and remove the 'serum' blade after the first load/application. Leave the 'urine' blade to complete the other 9 load/applications.

7. Following electrophoresis, (**SAS-1 Plus users:** remove the cover), remove the electrodes from the surface of the gel. (**SAS-3 users:** remove the alignment guide). Position the antiserum application template onto the gel surface. **NOTE:** The milled antisera channels should be aligned centrally over the printed box on the gel in which the samples are applied.
8. Apply 2 drops (or 50 μ L) of Protein Fixative (serum) or Total Antiserum (urine) into the hole of the SP lane and 2 drops (or 50 μ L) of the appropriate antiserum into the hole of the immunoglobulin lanes. Ensure that the fixative and antisera have completely filled the channels.
9. Incubate the gel. (**SAS-1 users:** incubate at 15...30°C).

10. Following incubation, place a blotter comb into the holes of the antiserum template. Allow 2 minutes for the excess antisera to be absorbed, then remove the blotter combs and the template from the gel surface.
11. **i) SAS-I users:** Remove the gel blocks using the Gel Block Remover and wash the gel in wash solution for 5 minutes.
ii) SAS-I Plus and SAS-3 users: Place a blotter D (smooth side down) onto the surface of the gel, leave for 10 seconds and remove.
12. **i) SAS-I users:** Place the gel on a blotter D agarose side up and place a blotter B (wetted in wash solution) onto the surface of the gel followed by two blotter X's. Press the gel using the Development Weight for 10 minutes.
ii) SAS-I Plus and SAS-3 users: Place a blotter D (smooth side down) onto the surface of the gel and replace the antiserum template to hold the blotter flat. Blot the gel.
13. **i) SAS-I users:** Remove the blotters and place the gel in wash solution for 4 minutes with gentle agitation.
ii) SAS-I Plus and SAS-3 users: Remove the blotter D.
14. **i) SAS-I users:** Remove the gel from the wash solution and place on a blotter D agarose side up. Place a blotter B (wetted in wash solution) onto the surface of the gel followed by a blotter D. Press the gel for 3 minutes.
ii) SAS-I Plus and SAS-3 users: Remove the gel blocks using the Gel Block Remover. Go to Step 16.
15. **i) SAS-I users:** Remove the blotters.
NOTE: Immediately after use, clean the antisera template with a mild, biocidal detergent. If possible, scrub the bottom of the template with a toothbrush or small test tube brush. Do not let the antisera dry on the template. Build up of antisera on the surface of the template will result in the formation of bubbles during the antisera application step. Dry the antisera template thoroughly. Water left in the holes will hinder the application of antisera for the next use. Store the template upside down to increase the air circulation and thus the drying potential of the applicator.
16. Attach the gel to the staining chamber holder.
17. Select the IFE test program on the staining unit and, following the prompts, Wash, Stain, Destain and Dry the gel.

a) SAS-2 (Auto-Stainer)

Step	Solution	Time (mm:ss)	Port	Temp (°C)
Dry	—	10:00		55
Wash	Wash solution	07:00	4	
Stain	Acid Violet Stain	03:00	5	
Destain	Destain solution	02:00	2	
Dry	—	05:00		65
Wash	Wash solution	03:00	4	
Wash	Wash solution	03:00	4	
Dry	—	05:00		65

b) SAS-4 (Auto-Stainer)

Step	Time (mm:ss)	Temp (°C)	Other
Wash	00:03		Recirculate ON
Wash	10:00		Recirculate ON

Stain	04:00	Recirculate ON
Destain	02:00	Recirculate ON
Destain	02:00	Recirculate ON
Dry	12:00	

63

c) **Manual**

Follow the sequence listed for the SAS-2 staining unit, using a staining bath for the Stain, Destain and Wash steps, and a Drying Oven with forced air at 60...70°C for the Dry steps.

- At the end of the staining cycle, remove the gel from the staining chamber. The gel is now ready for examination.

INTERPRETATION OF RESULTS

The majority of monoclonal proteins migrate in the cathodic, gamma region of the protein pattern, but due to their abnormal nature, they may migrate anywhere within the globulin region on protein electrophoresis. The monoclonal protein band on the immunofixation pattern will occupy the same position and shape as the abnormal band on the serum protein pattern. The abnormal protein is identified by the antiserum type it reacts with.

When low concentrations of abnormal protein are present, the abnormal band may appear as a band within the normal polyclonal immunoglobulin. A band can also be seen within a polyclonal background when there is a large polyclonal immunoglobulin presence also.

The publication 'Immunofixation for the Identification of Monoclonal Gammopathies' is available from Helena BioSciences on request.

Urine Immunofixation:

Type of Proteinuria

Normal urine
Glomerular

Bands Observed On Gel

Small albumin band
Albumin, alpha-1,
beta, gamma

Proteins Present

albumin
albumin, alpha-1
antitrypsin, transferrin,
gamma globulins
Retinol Binding Protein,
beta2-microglobulin,
alpha-2 microglobulin.
immunoglobulins,
free light chains

Tubular

alpha-1, alpha-2, beta

Overflow

gamma or variable

LIMITATIONS

1. Antigen Excess

Antigen excess will occur if there is not a slight antibody excess or antigen / antibody equivalence at the site of precipitation. Antigen excess in IFE is usually due to an excess of the immunoglobulin in the patient sample. Antigen excess is characterised by prozoning (unstained areas in the centre of the immunofixed protein band, with staining around the edges). A higher dilution of the sample should be used in this event to optimise the immunoglobulin concentration.

2. Non-Specific Precipitation in All Immunoglobulin Lanes

Occasionally a completed IFE plate exhibits a precipitate band in the same position in every pattern across the plate. This may result from:

- IgM monoclonal immunoglobulins.**

IgM monoclonal proteins can adhere to the gel matrix. A band will appear in all 5 antiserum lanes of the gel. However, where the band reacts with a specific antiserum for the heavy chain and light chain, there will be an increase in size and staining intensity of the band, allowing the immunoglobulin type to be identified. Additional dilution of the sample will improve the discrimination between the IgM-antibody reaction and the non-specific staining of precipitated IgM protein in other lanes, simplifying the diagnosis.

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Samples with high titres of Rheumatoid Factor or other immune complexes may show a preprecipitate band at the sample application point. Reducing the sample with DTT or β -2-mercaptoethanol can eliminate this non-specific reaction (Mix 190 μ L of diluted serum to 10 μ L of 1% (w/v) DTT in 0.85% saline solution or mix 100 μ L of serum with 10 μ L of a 1:10 dilution of β -2-mercaptoethanol in water. Perform the IFE as usual. Note: Always work in a fume hood when using β -2-mercaptoethanol).

c) Fibrinogen.

Fibrinogen, if present in the sample, will show as a discrete band in all lanes of the immunofixation pattern. Fibrinogen is present in plasma, and sometimes in the serum of patients on anticoagulant therapy.

3. Reaction With Kappa or Lambda Light Chain Antisera but No Reaction with IgG, IgA or IgM Heavy Chain Antisera.

Samples showing this pattern may either have a free light chain monoclonal gammopathy or they may have an IgD or IgE monoclonal protein. In this situation, the IFE should be repeated, substituting IgD and IgE antisera for two of the other heavy chain antisera. Failure to obtain a reaction with IgD or IgE antisera would be indicative of free light chain disease.

4. Band In Gamma Region Showing No Reactivity With IFE Antisera.

C Reactive Protein (CRP) may be detected in patients with acute inflammatory response¹²⁻¹³. CRP appears as a narrow band at the cathodic end of the serum protein pattern. Elevated Alpha1-Antitrypsin and Haptoglobin are supportive evidence for CRP. Patients with a CRP band will probably have an elevated level when assayed for CRP. A narrow band on the point of sample application can sometimes be seen which can be caused by chylomicrons in the serum or precipitated protein in samples which have been stored frozen.

5. Non-Reactivity With Kappa and Lambda Antisera

Occasionally a sample will have a reaction with a heavy chain antiserum but no light chain reaction is obvious. In this situation, the following need to be ruled out - a) Heavy chain disease, b) Very high concentrations of light chains, leading to antigen excess, c) Low concentrations of light chains, d) Atypical light chain molecule that does not react with the antiserum, e) Light Chains with 'hidden' light chain determinants (as sometimes seen with IgA and IgD). To obtain definitive results, testing may include a) A higher or lower dilution of the sample to optimise the antibody / antigen equivalence, b) Antisera from more than one manufacturer to aid in the identification of atypical immunoglobulins, and c) Treat the sample with β -2-mercaptoethanol to 'reveal' the light chains.

PERFORMANCE CHARACTERISTICS

A series of samples were tested and compared to another commercially available test kit - both kits showed equivalent results.

QUALITY CONTROL

Helena Biosciences IFE Control Kit (Cat. No. 9400) can be used to confirm the presence of monoclonal banding in all antisera lanes.

Helena Biosciences Kemtrol Abnormal Serum Control (Cat. No. 7025) can be diluted 1 in 100 and used as a positive control for Urine IFE.

BIBLIOGRAPHY

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980 ; May/June : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab.; 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Afonso, E., 'Quantitation Immunoelectrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
7. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox. Sang., 1969; 17 : 445-452.
8. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.
9. Johnson, A.M., 'Genetic Typing of Alpha(1)-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of Pi M.', J. Lab. Clin. Med., 1976; 87 : 152-163.
10. Cawley, L.P., Minard, B.J, Tourtellotte, W.W., Ma, B.I. and Chelle, C., 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
11. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
12. Jeppsson, J.O., Laurell, C.B. and Franzen, B., 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
13. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M., 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

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BIBLIOGRAFIA:

1. Fauchier, P. and Catalan, F. 'Interpretive Guide to Clinical Electrophoresis' Alfred Fournier Institute, Paris, France, 1988.
2. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M. 'Finding Clues to Disease in Urine' Diagnostic Medicine, 1980 ; May/June : 69-75.
3. Umbreit, A. and Wiedemann, G. 'Determination of Urinary Protein Fractions. A Comparison With Different Electrophoretic Methods and Quantitatively Determined Protein Concentrations' Clin. Chim. Acta., 2000; 297 : 163-172.
4. Wiedemann, G. and Umbreit, A. 'Determination of Urinary Protein Fractions by Different Electrophoretic Methods', Clin. Lab.; 1999, 45 : 257-262.
5. Wong, W.K., Wieringa, G.E., Stec, Z., Russell, J., Cooke, S., Keevil, B.G. and Lockhart, S. 'A Comparison of Three Procedures for the Detection of Bence-Jones Proteinuria' Ann. Clin. Biochem., 1997, 34 : 371-374.
6. Afonso, E., 'Quantitation Immunoelectrophoresis of Serum Proteins', Clin. Chim. Acta., 1964; 10 : 114-122.
7. Alper, C.A and Johnson, A.M., 'Immunofixation Electrophoresis: A Technique for the Study of Protein Polymorphism', Vox. Sang., 1969; 17 : 445-452.
8. Alper, C.A., 'Genetic Polymorphism of Complement Components as a Probe of Structure and Function', Progress in Immunology. First International Congress of Immunology. 1971 : 609-624, Academic Press, New York.

9. Johnson, A.M., 'Genetic Typing of Alpha(1)-Antitrypsin in Immunofixation Electrophoresis. Identification of Subtypes of Pi M.', J. Lab. Clin. Med., 1976; 87 : 152-163.
10. Cawley, L.P., Minard, B.J, Tourtellotte, W.W., Ma, B.I. and Chelle, C., 'Immunofixation Electrophoretic Technique Applied to Identification of Proteins in Serum and Cerebrospinal Fluid', Clin. Chem., 1976; 22 : 1262-1268.
11. Ritchie, R.F and Smith, R. 'Immunofixation III, Application to the Study of Monoclonal Proteins', Clin. Chem., 1976; 22 : 1982-1985.
12. Jeppsson, J.O., Laurell, C.B. and Franzen, B., 'Agarose Gel Electrophoresis', Clin. Chem., 1979; 25 (4) : 629-638.
13. Killingsworth, L.M., Cooney, S.K. and Tyllia, M.M., 'Protein Analysis', Diagnostic Medicine, 1980; Jan/Feb : 3-15.

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