### INTENDED PURPOSE

The Ristocetin Cofactor Assay kit is intended for the quantitative determination of an activity that reflects von Willebrand factor activity in plasma<sup>1</sup>.

The Ristocetin Cofactor Assay Kit measures the ability of a patients plasma to agglutinate formalin-fixed platelets in the presence of ristocetin. The rate of ristocetin-induced platelet agglutination is related to the concentration of von Willebrand factor.

### WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for In-Vitro diagnostic use only - **DO NOT INGEST**!. Wear gloves when handling all kit components. Refer to the product safety data sheets for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample. Some components of this kit contain sodium azide as preservative. To prevent the build up of potentially explosive metallic azides in metal plumbing, flush sinks with copious amounts of water and decontaminate regularly with sodium hydroxide solution.

### COMPOSITION

- Lyophilised Platelets (4 × 5ml)**

Contains formalin-fixed platelets in Tris-buffered saline.  
**Preparation:** Reconstitute each vial with 5.0ml of Tris-buffered saline. Swirl gently and allow to stand for 20 minutes.
- SARP Calibrator (2 × 1ml)**

Contains 1ml of lyophilised normal human plasma assayed for Ristocetin Cofactor activity.  
**Preparation:** Reconstitute each vial with exactly 1.0ml of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
- Tris-Buffered Saline (1 × 35ml)**

Contains 35ml of buffer containing 0.05M Tris and 0.15M sodium chloride with sodium azide as preservative. The buffer is ready for use as packaged.
- Ristocetin Cofactor Abnormal Control (2 × 0.5ml)**

Contains 0.5ml of lyophilised plasma from an individual with von Willebrand disease. The activity is printed on an insert sheet provided in the kit.  
**Preparation:** Reconstitute each vial with 0.5ml of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
- Ristocetin (2 × 1.5ml)**

Contains 10mg/ml Ristocetin lyophilised with stabilisers.  
**Preparation:** Reconstitute each vial with 1.5ml of purified water. Swirl gently and allow to stand for 10 minutes. Mix gently before use. Do not shake.
- Other kit components**

Each kit contains Instructions For Use, lot specific reference values insert and a report form.

### STORAGE AND SHELF-LIFE

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

- Lyophilised Platelets**

Once reconstituted, the reagent is stable for 1 day (24 hours) at 2...6°C or 30 days at -20°C.
- SARP Calibrator**

Once reconstituted, the reagent is stable for 2 hours at 2...6°C or 30 days at -20°C.
- Tris-Buffered Saline**

Opened bottle must be stored at 2...6°C.
- Ristocetin Cofactor Abnormal Control**

Once reconstituted, the reagent is stable for 8 hours at 2...6°C or 30 days at -20°C.
- Ristocetin**

Once reconstituted, the reagent is stable for 8 hours at 2...6°C or 30 days at -20°C.

### SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000 - 3000 xg for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

### STEP BY STEP PROCEDURE

Prepare all reagents as instructed with each pack. Set up platelet aggregometer as indicated in the instruction manual. For half volume settings see AggRAM step by step guide, document HL-2-1785P DOI: 2008/01 (1)

- Standard Curve Preparation**

Prepare a standard curve (in TBS) with suitable dilutions of SARP to give at least a 3 point curve between 20% and 90%. Mix without shaking.
- Patient Sample Preparation**

Prepare 1+1 and 1+3 dilutions of the patient plasma or control plasma in Tris-buffered saline. Mix without shaking.
- Aggregation Blank Preparation**

Pipette 250µl of platelet suspension and 250µl of Tris-buffered saline into an aggregometer cuvette with a stir-bar. This is used to set the 100% activity level.
- Testing**

- Set the 100% activity in each channel using the Aggregation Blank sample.
  - Pipette 400µl of platelet suspension into an aggregometer cuvette for each sample to be tested. Add a stir-bar to each cuvette.
  - Place the cuvette into the aggregometer and add 50µl of Ristocetin to the cuvette. Incubate at 37°C for exactly 3 minutes.
  - Add 50µl of standard, patient or control plasma dilution.
  - Monitor platelet aggregation for 5 minutes.

### INTERPRETATION OF RESULTS

Plot SARP dilution activity (X-axis) versus Aggregation Slope (Y-axis) on log-log graph paper. Draw a straight line of best fit. Interpolate patient slope values from the graph to give % activity. Correct these results for dilution and assayed activity of the SARP plasma.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

The Ristocetin Cofactor Abnormal Control supplied with the kit should be used undiluted and should give values within the range provided.

Helena BioSciences supply the following controls available for use with this product: REF 5301 SAC-1 REF 5302 SAC-2

### LIMITATIONS

The Ristocetin Cofactor Activity fails to reflect accurately von Willebrands' disease in several situations such as pregnancy<sup>1</sup>, infusion of commercial factor VIII concentrates<sup>3,4</sup> or administration of DDAVP5. In such instances, vWF R:Co may be corrected, yet the bleeding time remains prolonged. In addition, vWF R:Co levels may be normal in Type II B von Willebrands' disease even though the bleeding time is prolonged<sup>6</sup>.

### REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish it's own normal range. Expected values for Ristocetin Cofactor Activity in adults are 58-166% (0.58-1.66 units/ml). Each laboratory should determine its own normal values for the instrument / reagent system in routine use.

### PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena BioSciences or their representatives as a guideline. Each laboratory should establish it's own performance data.

The Helena BioSciences Ristocetin Cofactor assay is designed to give a linear calibration curve from 20 - 90%, without pre-dilution. Within run CV's are expected to be <10%.

### BIBLIOGRAPHY

- Zimmerman, T.S and Ruggeri, Z.M.'von Willebrands' Disease', Progress in Thrombosis and Haemostasis, Spaett (Ed), 1982, 6: 203-236.
- Ratnoff, O.D. and Bennett, B. 'Clues to Pathogenesis of Bleeding in von Willebrands' Disease', N. Eng. J. Med., 1973, 289 : 1182-1183.
- Blatt, P.M. et al. 'Antihemophilic factor Concentrate Therapy in von Willebrands' Disease' J. Am. Med. Assn., 1976, 236 : 2770-2772.
- Green, D. and Potter, E.V., 'Failure of AHF Concentrate to Control Bleeding in von Willebrands' Disease' Am. J. Med., 1976, 60 : 357-360.
- Mannucci, P.M. et al. 'Studies on the Prolonged Bleeding Time in von Willebrands' Disease', J. Lab. Clin. Med., 1976, 88 : 662-671.
- Ruggeri, Z.M. et al., 'Heightened Interaction Between Platelets and Factor VIII/von Willebrand Factor in a New Subtype of von Willebrands' Disease.', N. Eng. J. Med., 1980, 302 : 1047-1051.

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## BIBLIOGRAFIA

- Zimmerman, T.S and Ruggeri, Z.M.'von Willebrands' Disease', Progress in Thrombosis and Haemostasis, Spaett (Ed), 1982, 6: 203-236.
- Ratnoff, O.D. and Bennett, B. 'Clues to Pathogenesis of Bleeding in von Willebrands' Disease', N. Eng. J. Med., 1973, 289 : 1182-1183.
- Blatt, P.M. et al. 'Antihemophilic factor Concentrate Therapy in von Willebrands' Disease' J. Am. Med., Assn., 1976, 236 : 2770-2772.
- Green, D. and Potter, E.V., 'Failure of AHF Concentrate to Control Bleeding in von Willebrands' Disease' Am. J. Med., 1976, 60 : 357-360.
- Mannucci, P.M. et al. 'Studies on the Prolonged Bleeding Time in von Willebrands' Disease', J. Lab. Clin. Med., 1976, 88 : 662-671.
- Ruggeri, Z.M. et al., 'Heightened Interaction Between Platelets and Factor VIII/von Willebrand Factor in a New Subtype of von Willebrands' Disease.', N. Eng. J. Med., 1980, 302 : 1047-1051.

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- Zimmerman, T.S and Ruggeri, Z.M.'von Willebrands' Disease', Progress in Thrombosis and Haemostasis, Spaett (Ed), 1982, 6: 203-236.
- Ratnoff, O.D. and Bennett, B. 'Clues to Pathogenesis of Bleeding in von Willebrands' Disease', N. Eng. J. Med., 1973, 289 : 1182-1183.
- Blatt, P.M. et al. 'Antihemophilic factor Concentrate Therapy in von Willebrands' Disease' J. Am. Med. Assn., 1976, 236 : 2770-2772.
- Green, D. and Potter, E.V., 'Failure of AHF Concentrate to Control Bleeding in von Willebrands' Disease' Am. J. Med., 1976, 60 : 357-360.
- Mannucci, P.M. et al. 'Studies on the Prolonged Bleeding Time in von Willebrands' Disease', J. Lab. Clin. Med., 1976, 88 : 662-671.
- Ruggeri, Z.M. et al., 'Heightened Interaction Between Platelets and Factor VIII/von Willebrand Factor in a New Subtype of von Willebrands' Disease.', N. Eng. J. Med., 1980, 302 : 1047-1051.

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